

ACTIVESHOP[®]

Instructions for use

Saeyang Marathon 3 Champion milling machine



Importer:

ACTIVESHOP SP. Z O.O.
8b Graniczna St., building
DC2A 54-610 Wrocław,
Poland at
www.activeshop.com.pl
Rev. A/04/2024



Professional milling machine Marathon-3 Champion, is characterized by high quality workmanship, small dimensions and modern design. Clear and easy-to-read control panel, intuitive operation, creates an excellent opportunity for both experienced and novice users to work. The device is equipped with a very powerful head generating a torque of 2.9 NCm, so that already at 3000 revolutions it achieves maximum power. The vibration of the head has been completely eliminated, and the ergonomic shape and lightness will make working with it a pure pleasure.

The head has a TWIST-LOCK type cutter lock, which allows you to replace the cutter in one motion. Simply turn the handle (10) clockwise and replace the cutter, then in the opposite direction to lock the cutter.

Control of the milling machine is possible from the base and foot drive. An additional head holder can also be attached to the device.

BEFORE TURNING THE UNIT ON, CHECK THAT THE VOLTAGE LEVEL SWITCH IS SET TO 220V. SWITCHING ON A DEVICE SET TO 110V WILL BURN IT OUT.



<ol style="list-style-type: none"> 1. Speed control 2. On/off switch 3. Head plug socket 4. RPM switch 5. Control switch- hand/leg 6. Foot drive plug socket 7. Power cord 	<ol style="list-style-type: none"> 8. Voltage selection (Do not switch to 110V) 9. Micromotor 10. Handle 11. Cutter 12. Head stand 13. Cutter adjustment wrench 14. Foot drive 15. Head stand
---	---


Attention!

Unlocking the handle during operation will damage the machine. Changing cutters while the engine is running will damage the clutch. This damage is not subject to warranty. The cutter should be replaced in the open position of the handle (unlocked) and the cutter must be turned off at this time.

RECOMMENDED SPEED DEPENDING ON THE TYPE OF CUTTER

Cutter diameter Rotation max

5mm	41 000
7mm	30 000
10mm	20 000
13mm	16 000
16mm	10 000




NOTE!

The head should be cleaned and protected from dust.


Due to poor use and damage to the head due to dust, complaints will not be accepted.

Base



Technical data:

- power supply: 110-240 V/50-60Hz
- ergonomic housing
- clockwise and counterclockwise rotation
- on/off switch and speed control at base. Infinitely variable speed control from 0 to 35 thousand.
- twist-lock handle (replace cutter with one movement)
- cutter diameter 2.32 mm
- 2.9 Ncm torque - milling machine available in different colors



MODEL	DIMENSIONS IN mm			MASS
	Width	Length	Height	
Champion	118 mm	156 mm	78 mm	1300 g

HEAD	RPM	TORQUE	MOC	DIMENSIONS		MASS
				Length	diameter	
H20	0 - 30 000	2.9 Ncm	45 W	151	26	156 g
H200	0 - 30 000	2.9 Ncm	45 W	129	22	143 g
SH30N	0 - 30 000	2.9 Ncm	45 W	129	22	160 g
SH20N	0 - 30 000	2.7 Ncm	45 W	128	22	146 g
SH37LN	0 - 40 000	4.5 Ncm	45 W	152	28	210 g

H35SP1	0 - 35 000	3.2 Ncm	45 W	152	28	210
H37LSP	0 - 35 000	3.2 Ncm	45 W	154	30	215

Basic problems and solution:

The milling machine does not turn on.

- Check that it is plugged in - connect the milling machine to the power supply.
- Check that the milling machine is turned on - turn on the milling machine with the switch (No. 2 in the diagram).
- Check that the hand-leg switch (No. 5 in the diagram) is set to the "Hand" position.
- Check the head cable for damage - replace the head cable (available at Activeshop).
- Check that the power cable is not damaged - contact Activeshop service. - Check if the voltage change button (No. 8 on the diagram) is set to 220V. If it is not, and the milling machine is plugged in, the milling machine has been damaged - contact Activeshop service for repair against payment.

The head is overheating.

- Check that the Twist-Lock is properly closed.

Change the cutter after turning off the milling machine with the switch (No. 2 in the diagram). Turning the speed control potentiometer (No. 1 in the diagram) to the minimum does not turn off the milling machine, but only reduces its speed to the minimum (not visible to the naked eye). Changing the cutter on the switched-on milling machine at minimum speed can cause overheating of the head, and consequently damage (not covered by the warranty).

Attention!!! Do not soak the heads. Do not disinfect with excessive amounts of liquid. Water or disinfectant liquid getting inside the head, may damage it. Before possible disinfection, remove dust from the head. It is advisable to maintain the head with dedicated products.

SECURITY MEASURES

1. The device is a professional equipment for use in beauty parlors. Operation of the device only by qualified cosmetic personnel allows to avoid undesirable consequences during the treatment.
2. Do not disassemble the device or attempt to use it for purposes other than those indicated in this manual. Any service and repair work must be performed only by an authorized service technician (ACTIVESHOP Ltd. company).
3. Do not install or operate the device near liquids, wet places or with wet hands - so that no liquids get into the device.
4. Immediately remove the plug from the outlet and contact the authorized service provided by ACTIVESHOP Ltd. in the following cases:
 - a) liquid has entered the device, b) the device emits a strange odor, smoke or makes strange sounds, (c) the power cord was damaged, (d) the device was knocked over or its housing was damaged.
5. Do not put any objects on the cables. Do not install the unit where someone may step on the power cord.
6. For safety purposes, unplug the power cord from the outlet each time you finish working with the device.

7. Do not insert anything into the opening of the device. If a foreign object gets into the hole, remove the cord from the outlet and contact the distributor or manufacturer.
8. Do not place the device on an unstable cart, shelf or vehicle. If dropped, the device may be damaged.
9. In order to guarantee the established useful life of the device, it is recommended to keep an interval of one minute before pressing the power button again (between switching on and off and vice versa).
10. The power supply must comply with the specifications on the nameplate. Otherwise, there is a risk of incorrect operation of the device its damage or even burning.

NOTES:

1. In case of damage to the power cord or other malfunctions, contact an authorized service center - this will avoid danger.
2. This equipment may be used by children at least 8 years of age and by persons with impaired physical and mental capabilities and persons with inexperience and lack of familiarity with the equipment, if supervision or instruction is provided regarding the use of the equipment in a safe manner so that the associated risks are understood.
3. The device is not intended for play by children. Do not allow children to play with the device.
4. Do not leave the device unattended when children or persons without sufficient knowledge of how to use the device are near it.
5. This electrical appliance is not a toy. Improper use or contact with the electrical appliance may pose a health risk or lead to death, so use and store the appliance out of the reach of children.

DISPOSAL OF THE DEVICE



Dispose of used electrical appliances, fixtures and packaging to raw material reprocessing centers in accordance with environmental protection principles. Electrical appliances should not be disposed of in household waste containers. Follow the rules for discarding this type of equipment. In order to protect the environment in accordance with Directive 2012/19/EC on waste electrical and electronic equipment, used equipment

Electrical equipment should be collected separately, and taken to collection points for waste electronic and electrical equipment.

At the end of their useful life, as well as old appliances contain full-value r e c y c l a b l e components, so they should be handed over for recycling. Therefore, they should be handed over to appropriate collection points, as well as they can also be handed over to ACTIVESHOP Ltd.

STORAGE AND TRANSPORT

The device should be transported at a temperature in the range of (-) 20 °C to (+) 50 °C only in its original packaging. Store the device in a well-ventilated and relative humidity of no more than 75%. Operating temperature from (+) 10 °C to (+) 40 °C. During storage and transportation processes, the device must not be thrown or dropped, exposed to mechanical shocks or atmospheric agents.

